

FERSEMED

Heel Protector



It is used for heel protection decubitus ulcers (pressure sores) while the patient is reclining for a long term.



- It is used for heel protection decubitus ulcers (pressure sores) while the patient is reclining for a long term.
- Prevents the heel contact with sheets and prevents erosion of involuntary motion.
- Easy to use and simple.
- Wraps the heel comfortably and fixed by Velcro straps.
- Provides effective relief in patients.
- Can be washed at 30 °C and reuseable.

Fersemed Ordering Information

Model	Size* (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Fersemed	17 x 20	1	25	51 x 35 x 50	1750 gr	HPS 101