

PRODUCT INTRODUCTION CARD



Festermed Classic



Festermed Deluxe



Festermed Metal



Festermed QRB

Product Name:

Festermed

Definition:

Wrist / Ankle Fastening Bandage

Field of Use:

In turn, they have been produced for use in patients receiving treatment in Intensive Care units, patients who need to be controlled to prevent self-harm, patients on stretchers and cots who react very strongly.

In practice, they are used to fix the patient's hands and ankles to firm places such as cots, stretchers and so on

Function:

The tapes wrap the patient's hands and feet and fix them to the desired place (point) and restrict their movement.

Technical Specifications:

The products are designed in 4 different types for different places of use (see Model and dimensions to the title).

Almost all of the wide-width and robust-looking fastening tapes have velcro* locks on one end that serve to fasten the tape to the wrist after it is wrapped around the patient's wrist.

All of them also have connecting ribbons that allow the band to be fixed somewhereafter wrapping it around the wrists.

Except for the QRB model, all of the tapes are made up of 3 layers laminated to each other. The QRB model is produced by quilting method. The top layer consists of a chardonnayed ***(feathered) fabric. These hairs allow the velcro fixing particles to securely fix the tape.

There is a sponge layer in the middle part that provides softness to the tape.

The lowest layers of the bands that contact the skin are made of cotton combed fabric or terry fabric.

The products were sterilized by Gamma Ray for the first use and presented to the consumer, but;

Among them, the Classic model and the QRB model are disposable, but other models can be reused at 30 degrees Dec, in the washing machine, without washing and ironing by using household detergent.

(Velcro*: A special fabric -velcro fabric that is used instead of a snap or button,has micro-level plastic hooks on it and adheres to furry surfaces with these micro hooks).

(Laminating:** (Lamination) The process of gluing different materials together under heat or with the support of adhesive substances using special methods and apparatus and making them into one piece).

(Charding* :** Combing, combing the surfaces of cotton or synthetic fabrics with a special method the feathering process, while the chardonnay fabric is a fabric that has been made soft and velvety or terry-like by feathering one surface)

Model / Code / Measurements

It has four different models.

Festermed Classical is designed for patients in bedstead. It has two (2) different sizes for hands, foot,



For Hand: Code:ALH2001

it has a size of 25 x 5 cm and a thickness of 0.5 cm. The product has 2 connecting cords with a length of 75 cm.

For Foot: Code:ALH2002

it has a size of 40 x 5 cm and a thickness of 0.5 cm. it has 2 connecting cords with a length of 75 cm.,

Festermed Metal

This model is used for highly responsive patients on cots and stretchers..

Code: ALH 2011

It is 25 x 7 cm in size and 0.4 cm thick. It has a 150 cm long connection cord. This model is used to fix strong patients.

**Festermed Delüks**

This model is used for cases requiring strong fastening.

Code: ALH 2006

Measrement:

it measures 34 x 8 cm and is cm thick. The connecting cord is 25 cm, the length of the fixing cord, which is sewn to the end of the connecting cord and serves to fix it, is 32 cm, the thickness is 0.8 cm.

Festered QRB

is used in patients who are in a situation of blocking lifesaving treatment or further increasing tearing, pulling, scratching or peeling of the skin, are prone to self-harm or risk the integrity of the wound site. It fastens the patient to the bed with quick release buckels.



How to use;

In Festered Classic, Metal and Deluxe models, the bands are wrapped around the patient's wrist.

Velcro locking devices at the ends are attached to the bands and the bands are fixed to the wrist and finally the fixing cords are connected to the appropriate part of the stretcher or cot in order to restrict the movement of the wrist.

In the QRB model, the band is fixed to the wrist with a plastic buckle lock, and then the fixing cords are connected to the appropriate part of the stretcher or cot to restrict the movement of the wrist.

Caution:

That the bands are securely fixed to the wrists, do not interfere with blood circulation on the wrists,

It should be noted that the connecting cords are thoroughly connected to the fixing points.

Packing :

There is 1 product in each package, it is offered for sale in boxes of 25 pieces