


USAGE INSTRUCTIONS


FIXERBOARD

BABY WRIST HOLDER


MODEL : STANDART




Place the band on the wrist.



Place the band on the wrist.



Secure the connecting strips to the wrist with Velcro fasteners.



Check that the connecting strips do not obstruct blood circulation.



▷ **STERILE** Sterilized using gamma irradiation ▷ **Single Use**

▷ Do not use opened and damaged packs

Do not use opened and damaged packs

Do not use opened and damaged packs




MEDIKOKIM

USAGE INSTRUCTIONS


FIXERBOARD

BABY ANKLE HOLDER


MODEL : SMALL-MEDIUM-LARGE (LEG)




Place the band on the ankle.




Adjust according to the ankle.



Wrap the connecting strips of the band around the ankle.



Secure the connecting strips to the ankle with Velcro fasteners.



Check that the connecting strips do not obstruct blood circulation.



▷ **STERILE[®]** Sterilized using gamma irradiation ▷ **ⓧ** Single Use

▷ Do not use opened and damaged packs

©2017 MEDIKOKIM S.A. All rights reserved.

©2017 MEDIKOKIM S.A. All rights reserved. Do not use opened and damaged packs.