


USAGE INSTRUCTIONS

RESTMED


FREEDOM SPLINT

MODEL: STANDART/SOFT


REF NO : RSFRSW01- RSFRSW02- RSFRSW03- RSFRSW04- RSFRSW05-
RSFSSW01- RSFSSW02- RSFSSW03- RSFSSW04- RSFSSW05




Put the arm band on the patient's arm.



Pass the band sthrough the connecting buckles.



Secure the bands with Velcro fasteners. Check that the bands do not obstruct blood circulation.



▷ Non-steril

▷  Single Use

▷ Do not use opened and damaged packs