








Medikokim established in 1988 in Istanbul and started production process concerning its first product in 1990. In this period developed its product range and specifications.

The firm today; With its modern technology and expert staff, it continues to contribute to the health sector and the country's economy by producing class III, class I-s (Sterile) and class I-d (Sterile and without measuring function) disposable medical products.

In addition to Turkey, Medikokim products; It is used successfully in 30 countries in 5 continents, including Europe, Middle East, Asia, Africa, Australia, Central America and South America.











Index

Surgery & Operation	06	SAUGERMED Surgial Brain Pads
Products	09	CAPAD Drainage Cannula Pads
Intensive Care Products	12 13 14 16	NASOMED Nasal Cannula Fixation Bands TRAPAD Tracheal Cannula Pads TRAEGERMED Tracheal Cannula Holders STAENDERMED Endotracheal Cannula Holders
	18 18	PORTMED Multiple Catheter Holder PULSOMED Pulse Oximeter Probe Wrap
	19 20 21 22 23 25 26 26	FIXELMED Urine Bag Holder FIXERMED Male Urine Catheter Holder MITTMED Finger Control Mitts RESTMED Freedom Splint FESTERMED Hand & Foot Fixation Bands Fall Management Socks Fall Precaution Bracelets Essential Fall Management Kit
New Born Products	29 29 30 30 31	BINDERMED Baby Limb Holder NASOMED Baby Nasal Cannula Fixation Band FIXERBOARD Baby Limb Boards MAMASAFE Breast Protection Pad SCHÜTZERMED Phototherapy Eye Bands
Urology Products	35 36	CLAMPMED Urine Control Clamp CONDOFIX Condom Fixer
Wound Care Support Products	38 39	FERSEMED Heel Protector FERSEMED Rg Foot Elevator
Patient Safety Products	41 42 43 43 45 47 48 49	SAFEBELT Bed & Stretchers Safety Belts SAFEBELT Wheelchair Safety Belts SAFEBELT Patient Transfer Safety Belts SAFEBELT Optional Belts SAFEBELT Accessories SAFEBELT Wrist - Ankle Fixation Belts SAFEBELT Wheelchair Safety Belts SAFEBELT Patient Support Belts



Surgery & Operation Products








SAUGERMED

Surgical Brain Pads

For neurosurgical operations special produced haemostatic featured pads.



Enables safe use

Also, all the pads have an atraumatic featured Trace Rope with smooth touch. Pad locations can easily be traced and identified by these trace ropes. Pls. see figure 1

To avoid any distention Trace rope looped with a special technics*. Pls. see figure 2

All the pads consists of radio-pack blue strap which is X-ray detectable. Pls. see figure 3

To recognize the number of used pads during the operation it has a Special Counter Card 1 to 10. Pls. see figure 4

To prevent the pads from deformation it consists of a protective belt. Pls. see figure 5

Atraumatic

Produced from hydrofiber structured organic and inorganic materials in nonwoven form Atraumatic and soft.

Unfuzzy and unfeathered, does not leave loofah.

Sterile and dirtiness protected

Gamma sterilized.

Tvyek made of sterile pouch prevents unbounded tear and can only be open from its opening part.


It has an anticontamination bag in the sterile pouch.

There are 16 different dimensions of each pad made of 11 various raw material.

^{*} Patent rights of the special loop protected by the certification of the Turkish Patent Institute as Utility Model.









SAUGERMED

Surgical Brain Pads

Various options for different use of purpose 2 main featured and structured pad produced from 11 various raw material.

Viscose Pads





SAUGERMED Thinpad

- Absorbation capacity 7.1 times of its weight.
- Very thin, smooth and silky structured.
- Made of viscose and synthetic additive material.





SAUGERMED Plus

- Absorbation capacity 6 times of its weight .
- Thin, closely structured.
- Made of viscose and synthetic additive material.





SAUGERMED Thincot

- Absorbation capacity 8.5 times of its weight.
- Thin, smooth structured, laminated surface.
- Made of viscose and synthetic additive material.





SAUGERMED Laminat

- Absorbation capacity 8.3 times of its weight.
- Standard thickness, smooth structured, laminated surface.
- Made of viscose and synthetic additive material.





SAUGERMED Soft

- Absorbation capacity 7.3 times of its weight
- Standard thickness, closely structured, laminated surface.
- Made of viscose and synthetic additive material.





SAUGERMED Ultra Soft

- Absorbation capacity 8.5 times of its weight.
- Smooth , delicate structured.
- Made of viscose and synthetic additive material.






SAUGERMED Cotton Soft

- Absorbation capacity 10 times of its weight.
- Very smooth, delicate structured.
- Made of viscose and synthetic additive material.





Cotton Pads





SAUGERMED Pure Cotton

- Absorbation capacity 6.3 times of its weight.
- Very thin, very light structured.
- Made of 100% cotton raw materials





SAUGERMED Pure Cotton Thick

- Absorbation capacity 8.3 times of its weight
- Smooth, delicate structured.
- Made of 100% cotton raw materials





SAUGERMED Pure Cotton Medium

- Absorbation capacity 7.6 times of its weight
- Smooth, delicate structured.
- Made of 100% cotton raw materials





SAUGERMED Pure Cotton Soft

- Absorbation capacity 11.5 times of its weight
- Very soft, very absorbative, closely structured.
- Made of 100% cotton raw materials.

^{*}Please contact with our sales department for order informations of Saugermed.











FİXERBOARD-M
Wrist Support Bands



Intensive Care Products









Urine Bag Holder







Fall Management Socks





Essential Fall Management Kit



NASOMED

Nasal Cannula Fixation Bands

3 different model for adults, child and babies





NASOMED Adult

- Covers the head smoothly.
- Safe and easy to use with its velcro bands.
- Skin touch cotton surface absorbs the sweat
- In case of long term use does not cause allergy and irritation.
- Gamma sterilized.





NASOMED Child

- Consists of two parts supported both right and left ears.
- Safe and easy to use with its velcro bands.
- Can be used together with every size and model of cannula.
- Skin touch cotton surface absorbs the sweat.
- In case of long term use does not cause allergy and irritation.
- Gamma sterilized.





NASOMED Baby

- Practical and easy to use with its velcro bands.
- Keeps stable the cannula at the nose of the baby getting through the hoops of the nasal cannula.
- Can be used together with every size and model of cannula.
- Covers the head smoothly, skin touch cotton surface absorbs the sweat.
- In case of long term use does not cause allergy and irritation.
- Gamma sterilized.


Nasomed	Nasomed Ordering Information											
Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code						
Adult	65	1	25	11 x 20 x 25	355 gr.	NTH 5002						
Child	32	1	25	11 x 20 x 25	220 gr.	NTH 5001 TP						
Baby	43	1	25	11 x 20 x 25	245 gr.	NTH 5001						





TRAPAD

Tracheal Cannula Pads



3 different raw materials, 4 models from each raw material, Different sizes for each model for adults and children.

Practical

By its elastic cannula holes can be used with every type and model of cannula.

Atraumatic

Produced from smooth, hydrofiber structured of nonwoven or hydrophil foam material.

Offers alternative options

4 different models, each of them produced from various structured raw material. Each model has various dimensions for adult and children.

Sterile

Gamma sterilized.



Trapad Standart

Both sides absorbent.




Trapad Ultra

Double side pad, absorbent on both sides.




Trapad Blocker


One side absorbent the other side nonpermeable pad.




Straight slit Cannula hole at the top.



Without slit Cannula hole at the top.



Zigzag slit Cannula centered



Without slit Cannula centered.

(€₀₄₇₆

^{*}Please contact with our sales department for order informations of Trapad.



TRAEGERMED

Tracheal Cannula Holders

Two different strucured of holders: Sponge supported and rigid holders.

- Skin colour does not take attention.
- Skin touch cotton surface absorbs the sweat.
- In case of long term use does not cause allergy and irritation.
- Gamma sterilized.



Sponge Supported (Flexible) Fixation Holders

The sponge layer in the middle provides flexibility and softness to patient 2 different models.



Single piece model

Suitable for use in self.



Two pieces model

Better to use for the inpatient. Can be intervene to the cannula without moving.






Velor Supported Rigit Fixation Holders

Non flexsible and fixes canula in the first position thin structured fixation tapes 2 different models.



Single piece model

Proper for the use of the patients except inpatient.



Two pieces model

Better to use for the inpatient. Can be intervene to the cannula without moving.

Traegermed Sponge Supported Fixation Holders Ordering Information									
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code			
Single Piece Baby	23	1	25	11 x 20 x 25	200 gr.	TCHS 01			
Single Piece Child	45	1	25	11 x 20 x 25	230 gr.	TCHS 02			
Single Piece Adult	56.5	1	25	11 x 20 x 25	260 gr.	TCHS 03			
Two Pieces	37 + 6	1	25	11 x 20 x 25	290 gr.	TCHS TP			

Traegermed Velor Supported Fixation Holders Ordering Information										
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code				
Single Piece Baby	23	1	25	11 x 20 x 25	200 gr.	TCHR 01				
Single Piece Child	44.5	1	25	11 x 20 x 25	230 gr.	TCHR 02				
Single Piece Adult	56	1	25	11 x 20 x 25	260 gr.	TCHR 03				
Two Pieces	37 + 6	1	25	11 x 20 x 25	290 gr.	TCHR TP				

(€₀₄₇₆



STAENDERMED

Endotracheal Cannula Holders



CLP



STAENDERMED - CLP




Holder band fixes the cannula by gripping it with its locked clemp (CLP). There are 2 different models.

CLP - Single piece, velcro connection is on the back of the head model.


It grips and fixes the endotracheal tube with its clemp, it wraps around the neck and velcro s connects at the back of the head. There are 2 different sizes.

CLP TP - Two pieces, velcro connection is on the side of the head model.

It grips and fixes the endotracheal tube with its clemp, it wraps around the neck and velcro s connects at the side of the head. There is only one standard single size.









STAENDERMED - ST

Velcro Supported Holders

Holds the canula with the agency of velcro (V) lock, 3 different models

ST Single piece mid-supported

The middle portion of the tube wraps, 2 sizes.

ST TP Two pieces model

One size.




- Cotton surfaces touching the skin absorbs the sweat.
- Does not cause allergy and irritation in case of long term use.
- Cannula can be placed and removed without moving the patient's neck
- by the usage of single piece sidesupported and two pieces models.
- Three different types of Endotracheal Cannula Holders: Velcro supported, silicon supported and clemp supported.
- Gamma sterilized.





SCTP



STAENDERMED - SC

Silicone Supported Holder Bands

Holder band fixes the cannula by gripping it with its slicone surface (S). There are 3 different models.

SC <u>Single piece</u>, supported at middle part model.

The middle part wraps and grips the cannula. There are 2 sizes.

SCTP Two pieces model

There is only one standard single size.

Staendermed Ordering Information									
Model	Size (cm)		Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code		
ST	55	1	25	11 x 20 x 25	330 gr. ET	H ST 101	8698881471754		
01	65	1	25	11 x 20 x 25	345 gr. ET	'HM ST 102	8698881471761		
ST TP	50	1	25	11 x 20 x 25	330 gr. ET	TH ST TP	8698881472218		
SC	55	1	25	11 x 20 x 25	0	TH-SC 101	8698881478197		
	65	1	25	11 x 20 x 25	345 gr. ET	H-SC 102	8698881478333		
SC TP	43 + 30	1	25	11 x 20 x 25	330 gr. ET	TH-SC TP	8698881478326		
CLP	55	1	25	11 x 20 x 25	330 gr. ET	TH CLPM 101	8698881473734		
CLI	65	1	25	11 x 20 x 25	345 gr. ET	H CLPM 102	8698881473741		
CLP TP	50 + 30	1	25	11 x 20 x 25	330 gr. ET	TH CLP TP	8680508100932		

C€₀₄₇₆




PORTMED

Multiple Catheter Holder



- By its sewed catheter support band enables to fix the 3 different catheter at the same time.
- Easy to use, pratique and safe.
- Latex free.
- Can be fit to every arm size, does not hinder blood circulation.



- With the catheter pad dimensions of 15 x 8 cm, can be fix 3 different catheter at the same time.
- To distinguish the catheters easily it has 3 different colour of velcro binding.
- Does not cause allergy and irritation.
- Washable at 30 °C, reausable.

Portmed Ordering Information									
Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code			
Portmed	42 x 8.4	1	25	20 x 30 x 20	820 gr.	HP MED			

 ϵ

PULSOMED

Pulse Oximeter Probe Wrap



- Covers the probe safely and fix stable.
- Very simple to use.

- Can be used with every model and marked probe.
- Not sterile.

Pulsomed Ordering Information									
Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code			
Pulsomed	3 x 12	1	50	11 x 20 x 25	275 gr.	PPW 9001			





FIXELMED

Urine Bag Holder

Safely fix the urine bag holder to the leg and bears.





Practical

- A pair of fixation bands fits every size and model of leg urine bag.
- Safe and easy to use, can be washed at 30°C reausable.

Safe to use

- Elastic structure fits all leg size, covers the leg safely and keep stable.
- Wide part of the skin touch side which consists of silicone strap prevents the slide from the leg.

Comfort

- Special texture enables the skin breathe, reduces the sweat.
- Skin touch smooth feathery surface absorbs the sweat.
- Does not cause allergy and irritation.

Fixelmed O	Fixelmed Ordering Information									
Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code				
Fixelmed	60 x 2 / 60 x 24	1	10	11 x 20 x 25	320 gr.	UBH 1001				





FIXERMED

Male Urine Catheter Holder

Stables the urine catheter and prevents traumatic effects.





Practical

- Easy to use, pratique and safe.
- Fits all urine cannula types and models.
- Machine washable at 30 °C reusable.

Safe to use

- Latex free.
- Elastic structure fits all leg size, covers the leg safely and keep stable.

Comfort

- Special texture enables the skin breathe, reduces the sweat.
- Skin touch smooth feathery surface absorbs the sweat.
- Does not cause allergy and irritation.

Fixermed Ordering Information									
Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code			
Fixermed	60 x 5	1	25	11 x 20 x 25	680 gr.	CH 4001			

 ϵ



MITTMED

Finger Control Mitts

It is specifically designed for patients who try to remove catheters, attached cannulas or tend to hurt themselves.



MITTMED MITTFC 101 Fishnet

- It avoids the trauma with well rounded and soft wad part that placed hand palm.
- The upper part that consisting complete synthetic mesh restrict the movement of the fingers.
- It allows inserted use of IV cannula when get attached on the • Washable at 30 °C.
- hand because of its netting tissue.
- Easy and safe with its velcro bands.
- In case of long term using does not cause allergy and irritation.
- One model for both right and left



MITTMED MITTOTW Open

- It avoids the trauma with well rounded and soft wad part that placed hand palm.
- If necessary , can be interfered to the hand and fingers by opening the upper portion.
- Easy and safe with its velcro

bands.

- In case of long term using does not cause allergy and irritation.
- One model for both right and left arms.
- Washable at 30 °C.



MITTMED MITTCTW Closed

- It avoids the trauma with well rounded and soft wad part that placed hand palm.
- Easy and safe with its velcro bands.
- In case of long term using does not cause allergy and irritation.
- One model for both right and left arms.
- Washable at 30 °C.



MITTMED MITTSRVTW Therapy

- It helps to decrease patient stress as grabbing balls in the closed part of the mitt.
- Easy and safe with its velcro bands.
- In case of long term using does not cause allergy and irritation.
- One model for both right and left arms.
- Washable at 30 °C.

Mittmed Ordering Information										
Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code				
Fishnet	27 x 20	1	12	42 x 30 x 30	1100 gr.	MITTFC101				
Open	31 x 20	1	12	42 x 30 x 30	1600 gr.	MITTOTW				
Closed	32 x 20	1	12	42 x 30 x 30	1950 gr.	MITTCTW				
Therapy	29 x 22.5	1	15	42 x 30 x 30	850 gr.	MITTSRVTW				

CE



RESTMED

Freedom Splint

It is specifically designed for patients who try to remove catheters, attached cannulas or tend to hurt themselves. By preventing of the patient's elbow movement prevents hands touching through the head and body.



Standart Series

- 100% cotton lining; waterproof outer covering.
- In case of long term using does not cause allergy and irritation.
- It allows to using in place during X Rays.



- Easy and safe with its velcro bands and plastic clips.
- Clipped to the patient's sleeve to help prevent sliding.
- One model for both right and left arms.
- Washable at 30 °C.



Soft Series

- 100% cotton lining and soft, permeable to air, vinly mesh outer covering.
- In case of long term using does not cause allergy and irritation.



- It allows to using in place during X Rays.
- Easy and safe with its velcro bands and plastic clips.
- Clipped to the patient's sleeve to help prevent sliding.
- One model for both right and left arms.
- Washable at 30 °C.

Restmed Ordering Information										
Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Dimension (cm)	Box Weight	Order Code				
Standart XS	23 x 34	1	10	25 x 50 x35	1000 gr.	RSFRSW01				
Standart S	28 x 37	1	10	25 x 50 x35	1500 gr.	RSFRSW02				
Standart M	30 x 47	1	10	25 x 50 x35	2250 gr.	RSFRSW03				
Standart L	36 x 64	1	10	25 x 50 x35	4000 gr.	RSFRSW04				
Standart XL	49 x 69	1	10	25 x 50 x35	5000 gr.	RSFRSW05				
Soft XS	23 x 34	1	10	25 x 50 x35	1000 gr.	RSFSSW01				
Soft S	28 x 37	1	10	25 x 50 x35	1500 gr.	RSFSSW02				
Soft M	30 x 47	1	10	25 x 50 x35	2250 gr.	RSFSSW03				
Soft L	36 x 64	1	10	25 x 50 x35	4000 gr.	RSFSSW04				
Soft XL	49 x 69	1	10	25 x 50 x35	5000 gr.	RSFSSW05				

CE



FESTERMED

Hand & Foot Fixation Bands

Bands to fix the wrists of the patient in order to prevent the movement 5 different models for various purposes.





FESTERMED Metal

- For use of patients in bed or litter.
- With its perdurable connection ribbon through the metal hoop can be fixed to the needed place.
- Supported with its all round overlog seam extremely durable.
- Skin touch cotton surface absorbs the sweat.
- In case of long term use does not cause allergy and irritation.
- Can be used after washed 30 °C and making sterilization at autoclave.







FESTERMED QRB

- Patients assessed to be at risk of disrupting life-saving treatment or whose picking, pulling, scratching or peeling exacerbates a skin condition, causes self-injury, or compromises wound site integrity.
- Soft, and supported material.

- Initial hook and loop closure for additional security.
- Secures to the bed with quick-release buckles or ties.
- Adjustable size fits pediatric or adult limb.
- Wrist circumference 27 cm.

C€



FESTERMED Deluxe

- Especially for the patient who applies high strength.
- Skin touch cotton surface absorbs sweat.



- In case of long term use does not cause allergy and irritation.
- Washable at 30°C reusable.

C€





FESTERMED Classic

- \bullet For the use of inpatient both bed or litter.
- Skin touch cotton surface absorbs the sweat.
- In case of long term use does not cause allergy and irritation.
- Gamma sterilized.



Festermed Ordering Information										
Model	Size (cm)		Piece / Pouch		Pouch/ Box	Box Size (cm))	Box Weight	Order Code	
Deluxe	34	1	25	2	20 x 30 x 20	1215 gr.	/	ALH 2006	8698881473970	
QRB	25	1	25	2	0 x 30 x 20	1790 gr.	(QRB 2013	8682320930728	
Klasik El	25	1	25	2	20 x 30 x 20	850 gr.		ALH 2001	8698881471808	
Klasik Ayak	40	1	25	2	20 x 30 x 20	910 gr.	P	ALH 2002	8698881471815	
Klasik Bacak	65	1	25	4	20 x 30 x 20	1100 gr.	F	ALH 2003	8698881472577	
Metal El	25	1	25	2	20 x 30 x 20	1215 gr.	P	ALH 2011	8698881478425	



Fall Management Socks

Designed for fall risk patients.



- Patients requiring non-skid footwear to help prevent falls.
- Soft terry cloth.
- Skid resistant rubber tread on both sides.
- Available in 5 colors.


- 3 size: large, medium and small.
 - Small: up to 38 Medium: up to 45 Large: up to 47
- Latex free.
- 30 °C washable.

Fall Managem	Fall Management Socks Ordering Information								
Color	Size	Ordering Information							
	Small Socks Fits Up too Youth 38	MS35SR							
Red	Medium Sooks Fits Up too Adult 45	MS45MR							
	Large Socks Fits Up too Adult 47	MS47LR							
	Small Socks Fits Up too Youth 38	MS35SO							
Orange	Medium Socks Fits Up too Adult 45	MS45MO							
	Large Socks Fits Up too Adult 47	MS47LO							
	Small Socks Fits Up too Youth 38	MS35SP							
Purple	Medium Socks Fits Up too Adult 45	MS45MP							
	Large Socks Fits Up too Adult 47	MS47LP							
	Small Socks Fits Up too Youth 38	MS35SG							
Green	Medium Socks Fits Up too Adult 45	MS45MG							
	Large Socks Fits Up too Adult 47	MS47LG							
	Small Socks Fits Up too Youth 38	MS35SY							
Yellow	Medium Socks Fits Up too Adult 45	MS45MY							
	Large Socks Fits Up too Adult 47	MS47LY							



Fall Precaution Bracelets

Use for identify patients at risk of easily fall.



- Produced from the Polymer material compatible to the medical devices material terms and standards*.
- Antiallergic, antibacterial and antifungus hygienic structure that doesn't irritate the skin, not composed bacteria and fungal on it.
- Easy and rapid application.
- Safety lock with its poppered system.
- It has various sizes, length and width for different kind of purposes.
- *Not used cancerogenic and PVC/DOP featured materials harmful to human health.

Fall Prec	Fall Precaution Bracelets Ordering Information										
Model Dimensidon Color Size Piece/Box Box Size Box Weight Ordering											
DRHTB 1	Child	Yellow	2.5 x 24.2 cm	250	11 x 27 x 7 cm.	410 gr.	PSP - F01				
DRHTB 2	Adult	Yellow	3.2 x 30 cm	250	11 x 33 x 7 cm.	500 gr.	PSP - F02				



Essential Fall Management Kit*

- Easily identify patients at risk for falls.
- All components for high fall risk patient identification in one prepackaged kit.
- Takes the guesswork out of what products to use.

^{*}Please contact with our sales department for order informations of products.

FIXERBOARD-M

Baby Limb Boards



- Easy to use and pratique.
- Can be shaped to 3 dimension (downwards, upwards, right, left) by hand and can be stable at that shape.
- Both use for hand and leg wrist.
- Supports the wrists smootly, cotton structured skin touch surface absorbs the sweat.
- Does not cause allergy and irritation.
- Prevents the sprain of the wrist.
- With its connection ribbons and velcro lock enables to fix the wrist safely.
- Gamma sterilized.

Fixerboard Ordering Information										
Model	Size	Piece / Pouch	Pouch/ Box	Box Size (cm)		Order Code				
Small	6.5 x 14.5	1	25	11 x 20 x 25	288 gr.	LHC 04				
Medium	7.5 x 18.5	1	25	11 x 20 x 25	355 gr.	LHC 05				
Large	9 x 23	1	25	11 x 20 x 25	415 gr.	LHC 06				

C€₀₄₇₆



New Born Products













BINDERMED

Baby Limb Holder



- Easy to use and pratique.
- Both use for hand and leg wrist.
- Skin touch surface absorbs the sweat.



- Covers the wrist safely and stable with its connection ribbons.
- Does not irritate the skin and cause allergy.
- Gamma sterilized.

Bindermed Ordering Information										
Model	Model Size Piece / Pouch/ Box Box Order (cm) Pouch Box Size (cm) Weight Code									
Bindermed	3 x 11	1	25	11 x 20 x 25	263 gr.	BLH 3001				

C€₀₄₇₆

NASOMED

Baby Nasal Cannula Fixation Band



- Pratique and easy to use with its velcro bands.
- Can be used with every dimension and model of cannula.
- Covers the head smoothly, skin touch surface absorbs the sweat.
- Keeps the cannula safely at the babies nose getting through the hoops of the nasal cannula.
- Does not irritate the skin and cause allergy.
- Gamma sterilized.

Nasomed Ordering Information										
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)		Order Code				
Nasomed	38	1	25	11 x 20 x 25	245 gr.	NTH 5001				

(€₀₄₇₆



FIXERBOARD

Baby Limb Boards



- During the IV cannula application to the babies prevents the sprain of the wrist and support.
- Easy to use and pratique.
- Can be shaped to 3 dimension (downwards, upwards, right, left) by hand and can be stable at that shape.
- Both use for hand and leg wrist.
- Supports the wrists smootly, cotton structured skin touch surface absorbs the sweat.
- Does not cause allergy and irritation.
- Prevents the sprain of the wrist.
- With its connection ribbons and velcro lock enables to fix the wrist safely.
- Gamma sterilized.

Fixerboard Ordering Information										
Model	Size (cm)	Piece / Pouch/ Box Box Orde Pouch Box Size (cm) Weight Code								
Small	3 x 9	1	25	11 x 20 x 25	288 gr.	LHB 01				
Medium	4 x 10	1	25	11 x 20 x 25	355 gr.	LHB 02				
Large	4,4 x 12	1	25	11 x 20 x 25	415 gr.	LHB 03				

C€₀₄₇₆

FIXERBOARD-M

Baby Limb Boards



- During the IV cannula application to the babies prevents the sprain of the wrist and support.
- Easy to use and pratique.
- Can be shaped to 3 dimension (downwards, upwards, right, left) by hand and can be stable at that shape.
- Both use for hand and leg wrist.
- Supports the wrists smootly, cotton structured skin touch surface absorbs the sweat.
- Does not cause allergy and irritation.
- Prevents the sprain of the wrist.
- With its connection ribbons and velcro lock enables to fix the wrist safely.

Fixerboard Ordering Information											
		Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code					
Small	2.5 x 10.5	1	25	11 x 20 x 25	288 gr.	LHC 01					
Medium	3.5 x 12	1	25	11 x 20 x 25	355 gr.	LHC 02					
Large	4.5 x 13	1	25	11 x 20 x 25	415 gr.	LHC 03					

C€₀₄₇₆



SCHÜTZERMED

Phototherapy Eye Bands

- Prevents from the UV rays up to 450 Nm.
- Does not irritate the skin and cause allergy.
- Eye protective bands developed for the use of phototherapy treatment 4 different model.
- Gamma sterilized.





• Over head banded and supported model.



- Absorbtive cotton inner surface.
- 3 size: large, medium and small.



SCHÜTZERMED Sticker

• Sticker model adhesive to the temples.



• 3 size: large, medium and small.





SCHÜTZERMED Mask

• Model consists of head piece and protective mask.



• Adjustable for every size.



SCHÜTZERMED Yota

• Model consists of elastic band and protective eye cover.



- Net structure allows the UV rays to the head.
- 3 size: large, medium and small.

SchützermedOrdering Information											
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code					
Klasik Small	24.5	1	25	11 x 20 x 25	265 gr.	PEB 301					
Klasik Medium	30	1	25	11 x 20 x 25	298 gr.	PEB 302					
Klasik Large	32.5	1	25	11 x 20 x 25	320 gr.	PEB 303					
Sticker Small	11.5	1	25	11 x 20 x 25	240 gr.	PEB STC 301					
Sticker Medium	16.2	1	25	11 x 20 x 25	248 gr.	PEB STC 302					
Sticker Small	19	1	25	11 x 20 x 25	258 gr.	PEB STC 303					
Mask	19.2 + 25	1	25	11 x 20 x 25	325 gr.	PEB TPD					
Yota Small	23	1	25	11 x 20 x 25	218 gr.	PEB YT 301					
Yota Medium	27	1	25	11 x 20 x 25	240 gr.	PEB YT 302					
Yota Large	32	1	25	11 x 20 x 25	255 gr.	PEB YT 303					

(€₀₄₇₆

MAMASAFE

Breast Protection Pad

The product is placed onto breasts in order to protect nipples from trauma or any physical contact.





- The product consists of two different materials: Annular polyurethane foam and overlaying napped cotton woven fabric.
- The product is absorbent and soft.
- It absorbs milk leakage and sweat owing to its absorbent quality.
- It does not cause irritation on the skin and is antiallergic.

Mamas	Mamasafe Ordering Information											
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)		Order Code						
4 Set	Inner Dia: 5 h:3	1 Pair + Strapless	4 Pairs	26 x 10 x 11	172 gr.	BPP 104						
20 Set	Inner Dia:5 h:3	20 Pairs	20 Pairs	32 x 38 x 10	476 gr.	BPP 120						





Urology Products






CLAMPMED

Urine Control Clamp

The ideal solution is applied for 23 years in the incontinence problem!





Comfort

Body structure of Clampmed consists of semihard polymer of soft foam provides easiness and comfort during use.

With the balanced elastic bandage spreads its radial strength to the cellwall equally and reduces the clamp pressure to minimum level.

*Antiallergic featured of Clampmed is certified by the university report which is accredited by the Notify Body.

Safe

When applied to the penis special designed cushion at the bottom part pressured to the urine system to reduce the incontinence and does not effects the blood circulation.

Practical to use

In case of need easy to open without support and fix again.

Healthy and Hygienic

Sterilized for use, produced from antiallergic special materials so no irritation will be caused to the skin.* Also antibacterial and antifungal featured no bacteria and fungus will be composed**.

Confidant

Not occupy to much place unforseen easy to hide. Without thinking its presence very easy to use in daily life provides confidence to the patient.

Economic

Taking account of its using time, It is a more economical choice of product than absorbent

**Antibacterial and antifungus featured of Clampmed is certified by the laboratuar report which is accredited by the Notify Body.

Clampmed Ordering Information											
Model	Size	Piece / Pouch	Pouch/ Box	Box Size	Box Weight	Order Code					
	Normal	1	1	4 x 10 x 16 cm	31 gr.	UCC 01					
Clampmed	Large	1	1	4 x 10 x 16 cm	32 gr.	UCC 02					

CE



CONDOFIX

Condom Fixer



Practical

- Fixes the male external urine catheter safely and securely in place.
- Single size for adult and children.
- Non sterile and single use.

Condofix Ordering Information										
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code				
Condofix	3 x 12	5	10	15 x 25 x 21	325 gr.	CF 8001				





Wound Care Support Products







FERSEMED

Heel Protector

It is used for heel protection decubitus ulcers (pressure sores) while the patient is reclining for a long term.





- It is used for heel protection decubitus ulcers (pressure sores) while the patient is reclining for a long term.
- Prevents the heel contact with sheets and prevents erosion of involuntary motion.
- Easy to use and simple.
- Wraps the heel comfortably and fixed by Velcro straps.
- Provides effective relief in patients.
- Can be washed at 30 °C and reuseable.

Fersemed Ordering Information						
Model	Size* (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Fersemed	17 x 20	1	25	51 x 35 x 50	1750 gr	HPS 101

 ϵ



FERSEMED Rg

Foot Elevator

This foot elevation device is designed to reduce the risk of pressure ulcers (pressure sores) while the patient is in bed or reclining for a long term.





- It fits around the ankle to raise the heel off of any surface.
- It has a spiral-cut foam ring which easily fits over liner that provides foot elevation and comfort and can be fixed to the wrist with Velcro locks.
- Cylindrical body allows the patient's feet to the right to the left and back and forth.
- With its sponge structure provides softness and comfort.
- Cotton finished cover will not harm the skin.
- The product can be cleaned by hand or machine washed at 30 °C.

Fersemed Rg Ordering Information						
Model	Size* (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Fersemed Rg	R: 21 x h: 10	1	1	25 x 11 x 25 cm	304 gr.	HPC 101
	R: 21 x h:15	1	1	25 x 16 x 25 cm	344 gr.	HPC 102

 ϵ



Patient Safety Products



1 SAFEI

Bed & Stretchers Safety Belts



42

SAFEBELT



43

SAFEBELT





43

SAFEBELT



45

SAFEBELT



47

SAFEBELTBed Safety Belts Group



48

SAFEBELT
Wheelchair Safety Belts



49

SAFEBELT
Patient Support Belts



Patient Safety Belts

Steel Lock Serie



Safebelt patient safety belts used to control the movements of the patients with high safety risk (mental patients, alcohol addicts and patients at risk for escape) that might harm to their environments.

It has a wide product range to prevent and limited all undesired movements to respond all cases.

General Specifications

- By providing a secure fixation offers a maximum safety both for patients and around.
- Provides maximum freedom sensation to the patients where the patient movement is restricted as much as necessary.
- Application is easy and practical through effective magnetic switches.

- They are quite reliable with special steel pin locks.
- According to their usage areas, are divided into 3 main groups;
 - Used bad and stretcher
 - Used in chair
 - Used in patient transfer
- Products are made of 100% terylene and polyester material.
- 2 years warranty.
- Washable at 95 °C.





Bed & Stretchers Safety Belts

This group belts are used in order to fix patients securely to the bed and stretcher.



Waist Roll Belt

- Provides secure bedding for the patients.
- Provides freedom to turn to sides besides lying back and sitting positions.
- Meanwhile since patients are free to use their hands and legs, they can move comfortably for their social or personal needs (reading, eating etc.).





Shoulder Restraints

 Prevents shoulder and body movement and fixes the patients laid back. Used together with Waist Roll Belt or Multi - Function Extension Belt.



Wrist Restraints


- Securely fixes the wrist of the patient to the bed frame and restricts the movement certainly.
- Provides comfort for the patients due to its aesthetic and safe design.





Ankle Restraints

- Securely fixes the ankle of the patient to the bed frame and restricts the movement certainly.
- By using together with the extension, provides limited restriction for the patient's leg movement.
- Provides comfort for the patients due to its aesthetic and safe design.





Full Body Crisis Set

- Used to provide movement restriction for the patients during crisis.
- Consist of the items below;
 2 pcs. wrist restraints, 2 pcs. ankle restraints.
 Waist roll belt
 Shoulder restrains
 Multifunctional extension belts



Chair Safety Belts

This group belts are used in order to fix patients securely on the chair.





Wheelchair Waist Belts

- Secures the patient sitting on the wheelchair from slipping.
- If necessary can be used together with Pelvic Belt.





Wheelchair Shoulder Restraints

• Secures the patient sitting on the wheelchair from bending forward or backward also slipping.





Patient Transfer Safety Belts

Belt is used to restrict the patient's arm and leg movements while being taken from one place to another place. Belt is composed of three parts; including hand, foot and link belts.





Wrist Movement Restriction Belt

It is formed from two wrist belts and a belt wrapped around the waist and sewn onto the belt conveyor.

Wrist belts wrap securely and will not harm the skin ,conveyor belt ensure the hands at waist height by fixing the patient and prevents damage to itself and environment.





Ankle Restriction Belts

It is formed from conveyor belt and two wrist belts sewn onto the belt .

Wrist belts wrap the ankle comfortably without any pressure and in confidence that will not harm the skin.

The conveyor belt controls the patient's walking speed and stride length and prevents damage to the environment by restricting leg movement.





Binding Band for Wristbands and Ankle Bands

Binding band is used to join wristband set with ankle band set.



Support Belts

Supporting belts provides extra safety and comfort and used together with main safety belts.



Multi - Function Belt

Used together with wrist and ankle restraints to provide limited movement restriction.



Standard Extender

Used together with other main belts where longer ends needed.



Waist Extender

Provides extension for the waist belts for XXL size patients.



Pelvic Extender

Provides extra safety to the patients, used together with the waist roll belt and provides additional support in sitting position as well.

Safebelt Support Belts Ordering Information				
Model	Size*	Order Code		
W D. II D. I.	S-M (15-115)	SBB101		
Waist Roll Belt	L (15-145)	SBB102		
Wrist Restraints	Standart (14-26)	SBB201		
Wrist Restraints Set	Standart (14-26)	SBB202		
Ankle Restraints	Standart (16-32)	SBB301		
Ankle Restraints Set	Standart (16-32)	SBB302		
Ankle & Waist Restraints Set	Standart (14/16-26/32)	SBB605		
Shoulder Restraints	Standart (235)	SBB401		
Wheelchair Shoulder Restraints	Standart (216)	SBC 102		
Multi Function Belt	Standart (180)	SBS501		
Standard Extender	Standart (60)	SBS502		
Waist Extender	Standart (25.5)	SBS503		
Pelvic Belt	Standart (100)	SBS504		
Chandant Crisis Cat	S-M	SBB603		
Standart Crisis Set	L	SBB604		
Full Dady Crisis Cat	S-M	SBB601		
Full Body Crisis Set	L	SBB602		
Wheelchair Waist Belt	Standart (160)	SBC101		
Wrist Movement Restriction Belts		SBT201		
Ankle Restriction Belts	Standart (140)	SBT301		
Binding Band for Wrist Bands and Ankle Bands SBT501				



Accessories





Steel Lock

• Lock with magnetic control, special design, made of steel, can be unlocked by a special key.



Magnetic Key

• Easy to use, magnetically controlled key.



Bag

- Designed in 2 sizes can be able to carry complete set of crisis or small enough to carry single sets as separated.
- Existing of special cases which the locks and keys easily controlled.

Safebelt Accessoires Ordering Information				
Model	Size* (cm)	Order Code		
Steel Lock	Standart	SBA701		
Magnetic Key	Standart	SBA702		
Bag	Standart	SBA703		

(



Patient Fixation Belts

Buckled Soft Serie



Safebelt patient fixation belts used to control the movements in sensitive and precisely of the patients with secondary safety risk that might harm theirselves. It has a wide product range to prevent and limited all undesired movements to respond all cases.

General Specifications

- They provide maximum patient safety to itself and the environment by fixing in safe.
- They provide comfort and sense of freedom to the user even is restricted or controlled at the time of movement.
- According to their usage areas, are divided into 2 groups used in the chair and used in the bed.
- They are applied easily and practical because of their front fastening straps made of Velcro and shatter - resistant material made of POM buckled locks.
- Products are made of %100 terylene and polyester material.
- Can be washed at 95 °C.
- 2 years warranty.



Bed Safety Belts Group

This group belts are used in order to fix elderly patients in a safe and comfortably way to the bed.



Wrist Safety Belt



- Constants patient's wrist to the bed frame and restrict movement completely or partially.
- Provide comfort to the patients through its convenient model and aesthetic.
- Doesn't harm the skin cause to be made of soft material and provides comfort to the patient.



Ankle Safety Belt



- Constants patient's ankle to the bed frame and restrict movement completely or partially.
- Provide comfort to the patients through its convenient model and aesthetic.
- Doesn't harm the skin cause to be made of soft material and provides comfort to the patient.



Wheelchair Safety Belts

This group belts are used in order to fix elderly patients securely and comfortably to the wheelchair.





Wheelchair Safety Belt

- Prevents sliding and enables leaning forward for the patient sitting in the wheelchair.
- Provides support between a patient's leg with soft filled portion in the front.

Safebelt Bed Safety Belts Ordering Information				
Model	Size* (cm)	Order Code		
Wrist Fixation Belt	Buckled (45 cm)	YBM 201		
Ankle Fixation Belt	Buckled (65 cm)	YBM 301		
Wheelchair Fixation Belt	Buckled (164 -146 cm)	YBS 102		

 ϵ



Patient Support Belts

Safebelt Patient Support Belts are designed for patients who need support and at risk of falling.







Walking Belts

- Walking belt helps to support the patients who have difficulty to walk alone.
- Soft and durable belt has got six handles as being sewn on the vertical and horizontal ergonomic grip straps are available.
- By assistance of those handles on the belt wrapped above the waist helps prevent caregiver back injuries and assists in the safe transfer, stumble and fall risk is avoided.
- The belt is easily to grasp the patient's waist by its enough width.
- By courtesy of plastic buckle can be adjustable according to each body, opening and closing is simple, fast and secure.
- It can be cleaned by hand or machine washed at 95 °C.





Lift Assist Belt

- The belt used to lift up the patient from bed or wheelchair.
- Ergonomic hand grips are mounted at both ends of the
- The belt is wrapped around the patient's waist and given support to try to pull up, go back (fall risk) is blocked.
- It was made of durable soft, special flexible neoprene material.
- It can be easily carried by each officer because of lightness and takes up little space.
- If necessary, hand or machine washable at 30 °C.

Safebelt Patient Support Belts Ordering Information				
Model	Size* (cm)	Order Code		
Walking Support Belt	Standart (140 cm)	SBS 101		
Patient Lift Assist Belt	Standart (46 cm)	SBL 201		





Medikokim Tıbbi Kimyevi Malzeme San. Tic. AŞ İkitelli O.S.B. Demirciler Sanayi Sitesi B/7 Blok No:157 Başakşehir PK 34306 İstanbul / TÜRKİYE

Phone: +90 212 549 57 78 • Fax: +90 212 549 57 79

E-mail: info@medikokim.com

Rev.: 03