

TRAPAD

Tracheal Cannula Pads

CE 1783

3 different raw materials,
4 models from each raw material,
Different sizes for each model for adults and children.



Practical

By its elastic cannula holes can be used with every type and model of cannula.

Atraumatic

Produced from smooth, hydrofiber structured of nonwoven or hydrophil foam material.

Offers alternative options

4 different models, each of them produced from various structured raw material. Each model has various dimensions for adult and children.


Sterile

Gamma sterilized.



Trapad Standart

Both sides absorbent.




Trapad Ultra

Double side pad,
absorbent on both sides.




Trapad Blocker


One side absorbent
the other side
nonpermeable pad.




Straight slit Cannula
hole at the top.



Without slit Cannula
hole at the top.



Zigzag slit
Cannula centered.



Without slit
Cannula centered.

*Please contact with our sales department for order informations of Trapad.

