

## PRODUCT TRAINING CARD



### **Product name**

Condofix

### **Definition:**

Fixing Wrap

### **Field of Use:**

It is manufactured for Urology and Intensive care service. It is designed for male patients .

It is used for safely fixing the probe with a condom to the penis.

### **Function:**

After the condom probe is attached to the penis, the band is wrapped around the upper edge of the condom probe and prevents the condom probe from sliding down with the pressure it makes when it is tightened.

### **Technical Specifications:**

It has the appearance of a ribbon. At one end of the band there is a fixing accessory made of velcro\*, which allows the band to fix it when it gets entangled in the condom (penis).

The tape consists of 2 layers laminated\*\* to each other.

The top layer is composed of charonnayed \*\*\*(feathered) fabric. These hairs allow the velcro fastening apparatus to secure the tape safely. There is a sponge layer at the bottom that provides softness to the tape.

The product is disposable and non-sterile.

**Velcro\*:** A special fabric -velcro fabric that is used instead of a snap or button,has micro-level plastic hooks on it and adheres to furry surfaces with these micro hooks).

**Laminating\*\*:** (Lamination ) The process of gluing different materials together under heat or with the support of adhesive substances using special methods and apparatus and making them into one piece).

**Chardoning\*\*\*** : Combing, combing the surfaces of cotton or synthetic fabrics with a special method  
the feathering process, while the chardonnayed fabric is a fabric that has been made soft and velvety or terry-like by feathering one surface)

**Model**

The product is a single model.

**Code;**

CW 8001

**Measurements**

it is 12 cm long by 3 cm wide by 0.4 cm thick.

**How to use;**

The tape is wrapped on a condom attached to the body of the penis in an un-erect state, with a sweet firmness and secured with velcro.

**Caution:**

During the fixing process, care should be taken to ensure that the tape is tight enough to prevent blood circulation.

**Packing;**

There are 5 products in a bag, they are offered for sale in boxes of 10 pieces.